




GE Healthcare



Innova IGS 6

IGS 630 Configuration
Biplane Image Guided System
For Interventional Neuroradiology
and Combo use



Innova IGS 6 with AutoRight™

Advanced visualization for your evolving practice



PERFORM FAST INTERVENTIONS

Ergonomics, advanced anti-collision system and unique offset C-arm allow rapid Groin-to-Reperfusion time during mechanical thrombectomies.



HIGH DEFINITION 3D IMAGING FOR NEURO INTERVENTIONS

Understand complex vasculature, rule out hemorrhage and identify large vessel occlusions using high resolution CBCT and advanced 3D visualization tools.



GET THE RIGHT IMAGE AT THE RIGHT DOSE AUTOMATICALLY

AutoRight™^{1a} is the 1st AI-based interventional image chain™ in the industry, trained on 6,000+ datasets for automatic IQ-Dose optimization. Benefit from high DQE detectors, tailored to your clinical specialty.



IMPROVE OUTCOMES WITH AUGMENTED REALITY GUIDANCE

Perform any procedures while significantly reducing radiation dose and contrast media with intuitive ASSIST solutions.² Maximize clinical tools adoption with GE Clinical consultants.



GE HEALTHCARE SERVICE, DELIVERING PEACE OF MIND

Robust and reliable system with up to 99% guaranteed uptime³ and up to 38% reduction in unplanned downtime.⁴



Built on reliable foundations, the Innova IGS 6* interventional biplane X-ray system delivers innovative 2D and 3D imaging solutions for the head, the neck and spine.

With advanced applications and the capacity to customize, it helps you master the most complex cases, while a comprehensive set of versatile features can enable your lab to expand its mix of procedures to other interventional radiology and cardiovascular applications

So you have the right tools. For each patient. Every day.



Simplify your workflows

Exclusive Gantry Ergonomics

- ✔ The off-set position of the C-arm provides maximum positioning flexibility and patient access for coverage from groin to head, without the need to move the L-arm to the side.
- ✔ Thanks to the exclusive design of the C-arm, you get unrestricted access to the patient's head for anesthesia and nursing even when lateral plane is in place.
- ✔ The C-arm design allows to perform 3D rotational acquisitions of the spine and torso as well as patients' head while intubated.
- ✔ A detector size optimized for neuro interventions. With its 31 cm (12.2 in) square detector on each plane, there is no compromise on image quality between the planes. It maximizes the anatomical coverage while allowing to reach the desired working angles.



The Innova IGS 6* features dedicated biplane gantry design, customization and automation capabilities to accelerate procedures. **So you can focus on the patient.**

Like all GE systems, Innova IGS 6* is designed from the ground up using the tried and trusted GE imaging chain. It is optimized to provide the image clarity you need while helping you keep dose as low as possible.

Empowered with the multi-modality Advantage Workstation, plan, guide and assess procedures with confidence using dedicated advanced applications.

User interface meant for simplicity

- ✔ Control your system and images with integrated, intuitive and robust table-side controls
- ✔ With simple menus, the Central touchscreen lets you control the system, modify imaging parameters, configure your large display monitor and manipulate advanced applications – all at table-side.
- ✔ The full-color, 148 cm (58 in) medical grade large display lets you view multiple images from up to 16 sources, using up to 120 customized layouts.† Zoom in comfortably without loss of detail or pixilation to get the clinical focus you need.
- ✔ Save up to 63 positions (7 positions for each 9 users) on-the-fly from table side with the press of a button. The Auto-positioner allows you to reach gantry positions with the click of a button, from table-side or from the control room.

With one hand on a single console, drive both planes simultaneously or independently, both detectors, and the table motions

Control dose strategy at table side

Set up the system with the touch of a button thanks to user specified settings and presets protocols

Operate our advanced applications, third party connections and display options directly with a fingertip

Get more done. At table side. With the touch of a button



Table optimized for endovascular work

- ✔ Innova IGS 6* comes with the Omega V table, which allows +/- 180° table rotation around its vertical axis.
- ✔ The 333 cm (131") long table top accommodates catheters and provides additional area for placement of sterile items.
- ✔ It allows a total load of 304 kg (780 lbs.) to accommodate heavy patients up to 204kg (450 lbs.).
- ✔ With table panning, imaging coverage is as long as 195 cm (76") allowing to image from head to foot for the vast majority of patients.

Accelerated gantry motions and faster positioning

- ✔ Sensors on the frontal plane position the frontal detector automatically to avoid collisions with patient and the other plane.
- ✔ If CPR or surgery must be performed in emergency, total patient access can be provided at any time.
- ✔ Thanks to a gantry L-arm rotating at +/- 100° around its vertical axis, patients can be accessed easily from head to toe, from their right or left sides.
- ✔ Off-isocenter imaging - Each plane serves a purpose. The lateral plane can be iso-centered or off-iso-centered to allow visualizing vessels from two different centering positions with one injection.
- ✔ The combination of movements of the C-arm and the L-arm permits +/- 55° cranial and caudal angulations.
- ✔ Off-isocenter acquisition helps visualize both the origin of the vessel (lateral) and the distal area (frontal/AP) with just one injection and one acquisition.

Support that never sleeps

Your partner to help optimize patient care

1 Expert Service Delivery

- ✔ Customized and flexible offerings with service interventions guaranteed within 30 min remotely & within 4 hours on site, to achieve up to 99% uptime.³
- ✔ Dedicated and highly trained local experts, with 3830 modality-trained field engineers, 130 online engineers, and 160 customer agents in 17 centres in Europe.⁵
- ✔ Remote diagnostic solutions 24/7 (InSite, iInq) leading to 27% remote fix and 81% issues fixed in less than one visit on site.⁶

3 Education and Training

- ✔ Optimize equipment performance with customizable clinical application training with on-site, remote, and online options.
- ✔ Over 350 accredited continuing programmes and over 100 application specialists in Europe.⁵
- ✔ Connect to the GE Cares Community and learn how to capitalize on online resources to facilitate trainings and increase your professional skills.

2 OnWatch⁷ Convert unplanned to planned

We created OnWatch to maximize your efficiency by helping to ensure that your angiography system is operating when you need it to. OnWatch service measures key parameters from your equipment. It looks ahead to help limit disruption from unplanned downtime, creating a less stressful experience for you, your staff and your patients. This visionary technology drives progress in patient care, enhances efficiency and can help minimize the costs associated with downtime.



Join us now on GECARES.COM

Connect with other healthcare professionals and grow your network. Interact with key opinion leaders and view their publications.

Share your experience, publish content and stay up to date with the latest clinical trends shared by your peers.

Learn new techniques and increase your skills in your daily practice. Access online trainings, educational contents, clinical webinars built by experts for experts.

AutoRight™

AutoRight™, Automated Image Acquisition


6P with embedded neural networks (NNs)

- ✓ Up to 6 Parameters device control
- ✓ Automatic optimisation of acquisition parameters in real time, based on the imaged anatomy (6P)

AutoRight helps remove the burden of manual adjustment, lets you focus attention on patients.

Helps focus your attention on patient

- ✓ AutoRight is the **1st AI-based**, interventional image chain in the industry, trained on **6,000+** datasets
- ✓ Automatic adjustment of up to **6 parameters** in real time to optimise image quality and dose
- ✓ **2/3** of hardware and software renewed in the image chain




Powered By
Edison

AutoRight™, Intelligence inside

Providing you with the right image at the right dose, automatically

Automated image acquisition with dose and IQ optimisation in real time, whatever the patient anatomy or C-Arm angulation

Smart image processing to retain and present only the useful part of the information embedded in the acquired images.




Clear image display with high luminance DICOM monitor that is digitally linked to the image processor.

myIQ: you decide what's best for you!8

Image quality is not one-size-fits all for physicians.

On the contrary, preferred « looks » vary greatly from one clinician to another: myIQ lets you decide what's best for you! Configure your preferred looks with your local GE application support and get the best image quality for you. Where no two people are the same, why would your image system need to be?



Genuine look with limited digital processing. GE historical DSA look for devotees.

Respectful of X-ray characteristics, maintain small detail visibility, while reducing background noise.

Filters out most of the noise, ideal to image proximal and large vessel anatomies.

Sharp small details bright and contrasted.



Create the image You need

Digital Subtraction Angiography

Enhance visibility of vessels

The capability to image fine anatomical details is at the core of interventional neuroradiology. The Innova IGS 6 allows to reveal vascular structures through segmented DSA. Auto pixel-shifting allows to minimize impact of patient motions on image quality. Frame rate can go up to 7.5 fps for both planes simultaneously.

Improved de-noising and sharpening algorithms constantly enhance image quality in real time and after post-processing, without any compromise on the imaging plane.




Digital Subtraction Angiography

Single Shot

Enhance visibility of intracranial devices


Single Shot customizes the GE exclusive Digital Range Management to highlight the devices on flattened background.

Optimized for spatial resolution this acquisition helps you see the devices during their deployment with one foot pedal press.



Single Shot

Roadmapping applications can help you advance wires, catheters and devices through neurovascular structures, correcting for patient movement and performing on-the-fly adjustment of vessel transparency and pixel-shifting, all while minimizing contrast injections and X-ray radiation. The Innova IGS6 features 2 dimensional roadmapping modes with fluoroscopy and blended roadmap, or 3 dimensional with Vessel ASSIST¹⁴. Subtracted fluoroscopy is ideal for visualization of liquid embolic material, by providing full subtraction upon each new fluoroscopy pedal press.



Blended Roadmap

Navigate clearly in complex anatomy

This fast, easy, and potentially dose - and contrast - saving roadmapping application superimposes any Digital Subtracted Angiography image with 2D fluoroscopy, helping you view the advancement of guidewires and devices through vessels. Realtime pixel shifting quickly corrects for patient movement. Vessel transparency and image landscaping can be separately adjusted on the fly during fluoroscopy, all from the Central touchscreen. Blended roadmap also automatically adjusts to change of field of view.




Blended Roadmap

AngioViz

Visualize vascular flow

AngioViz facilitates understanding of vascular flow by displaying a DSA series in a single color-coded image of peak opacification, time to peak and combinations of those. AngioViz automatically synchronizes different DSA series for flow comparison of pre- and post-interventional runs.



AngioViz

Advanced

Applications

Empowered with the multi-modality Advantage Workstation, perform interventions with confidence using efficient, integrated and easy to use advanced applications.

Advantage Workstation Multimodality visualisation, analysis & navigation


The Advantage workstation (AW) VolumeShare has been created as a multimodality review workstation. By integrating datasets from interventional systems, as well as CT, MRI and PET, it provides remarkable convenience to compare and fuse image information. Now, interventional radiologists can perform processing, integration and image overlay at a single workstation with one user interface, even remotely with the AW server option.



Integrated Registration† Multi-modality image management

Integrated Registration²⁰ lets you fuse and register two volumetric acquisitions from the same or different acquisition modalities. With it, you can easily compare 3D anatomical images from CT, MR with PET, SPECT, and X-ray angiography⁹ for a more comprehensive analysis, and overlay them with live fluoroscopy.

Quickly compare 3D anatomical images from CT, MR PET, SPECT, and 3DCT HD¹⁰ by fusing two volumetric acquisitions and up to 20 different volumes.




Unprecedented 3D images in the neuro suite

3DCT HD¹⁰ Next-generation 3D imaging

Interventional procedures are becoming increasingly complex and clear 3D visualisation of soft tissue, anatomical structure and small interventional devices is crucial. CT-like imaging made available for interventional procedures has the potential to address even the most formidable challenges. 3DCT HD provides exceptional fine image details and quality on cross-sections to help you clearly visualise soft tissues, bones and small devices during your interventional procedures.

	3DCT HD
Spin duration	5, 7, 13 sec.
Frame rate	50 fps
Reconstructed 3D model resolution	512x512x512 256x256x256

After intra-op rotational acquisition, the Volume Rendered 3D and cross-sectional images are automatically reconstructed and displayed for planning and immediate assessment, with full 3D review from within the sterile field using the in-room wireless AW mouse.¹¹



Metal Artefact Reduction¹²

3DCT HD MAR¹² reduces streak artifacts induced by the presence of small metallic devices such as coils or clips within the 3D field of view.

Subtracted 3D¹³

Subtracted 3D enhances 3DCT HD by producing subtracted 3D vascular images from automated sequential mask and contrast spin acquisitions. Clinicians may use Subtracted 3D to quickly visualize vessels without the need to remove surrounding bone, tissue, and implanted devices. Interventional devices such as coils, stents, glue, and clips, as well as calcified plaque, are visible on the mask image and can be fused onto the subtracted image.




Vessels

Higher spatial resolution and automatic streak artifact reduction to image small vessels



Devices

Higher spatial resolution & Automatic reduction of metallic artifacts to visualize devices



Soft tissue

Scatter reduction for flatter images and higher tissue differentiation

Ischemic stroke thrombectomy

GE helps stroke teams achieve fast reperfusion of acute ischemic stroke patients and efficient secondary prevention supporting them with optimization of operational processes, innovative imaging, digital and education solutions



Vessel ASSIST

Vessel ASSIST¹⁴ delivers accurate and easy anatomy segmentation and accurate vessel quantification from 3D volumes.

It also offers advanced 3D roadmap and augmented reality capabilities to guide catheter with confidence, whether you wish to perform a stenting of an intracranial aneurysm with an intraluminal or intrasaccular flow diverter, peripheral chronic total occlusion recanalization, a uterine fibroid embolization, a tumor embolization prior to laparoscopic partial nephrectomy or a prostate artery embolization among others.



Vessel ASSIST enables:



With 3DCTHD and Virtual Dilation visualize the relationship between the vessels and the most complex devices such as intrasaccular or intraluminal stents, without the need to dilute your contrast injections.

WITH ASSIST FOR
VENOUS STENTING,
USERS ACHIEVED

a reduction of
the occurrence of delayed
flow in the cortical veins¹⁵



WITH ASSIST
USERS EXPERIENCED

an increase of operator's
confidence to identify
critical vascular structure
& device position¹⁵

VESSEL ASSIST ENABLES

robust and high quality
assessment
of the stent-to-vessels
relationship

Aneurysm embolization with Flow Diverter stent



PLAN

Segment and size the anatomy of interest in CBCT¹. Plan the stent centerline and landmarks and overlay them on the 3D model and prepare 3D landmarks to help guide the deployment.

GUIDE


Export segmented anatomy, centerline and landmarks and overlay them on live fluoroscopy on the frontal or lateral plane for augmented 3D fusion guidance.

ASSESS

With 3DCTHD² and Virtual Dilation³ visualize the relationship between your device and vessel with robust image quality and without changing your contrast injection protocols.

1. Requires VesselIQ Xpress, part of Vessel ASSIST.
2. 3DCT HD is an option sold separately. Includes 3DxR. Requires: AW Workstation and Volume Viewer. Available only on Innova IGS 5, Innova IGS 6, Discovery IGS 7 and Discovery IGS 7 OR.
3. Virtual Dilation is a GE created customized protocol which requires Volume Viewer and Integrated Registration. These applications are sold separately.

MR-guided Cerebral Venous Sinus Stenting



PLAN

Segment and size the anatomy of interest with 2 clicks in MR venogram.⁴ Plan the stent positioning on the 3D model and prepare 3D landmarks to help guide the deployment.

GUIDE

Import segmented anatomy, vessel contours and landmarks and overlay them on live fluoroscopy on the frontal or lateral plane for augmented-MR biplane fusion guidance.

ASSESS

Confirm stent positioning and venous blood flow in biplane DSA runs.

4. Require AW workstation with Volume Viewer and Volume Viewer Innova.

Embo ASSIST With Virtual Injection

Simulate your embolisation procedures with Virtual Injection

Embo ASSIST[®] with Virtual Injection, powered by Edison, is a 3D Visualization software solution designed to help clinicians simulate injections dynamically and thus perform embolization procedures with confidence. It offers simulation capabilities for complex procedures such as arterio-venous malformation, distal aneurysm or tumor embolization.



Embo ASSIST enables:



using 3D overlay of selected vessels




of vasculature from a CBCT acquisition



the destination of potential injections

Arterio-Venous Malformation embolization



PLAN

With 3DCT HD¹ and Split Phase², visualize arterial and venous phases details. Developed to automatically track vessel pathways to anticipate navigation challenges. Designed to allow you dynamically test embolization strategies with Virtual Injection and define optimal catheter positions.


GUIDE

Export automatically segmented vessels and landmarks of potential embolization points for augmented-3D fusion guidance. Overlay them on live fluoroscopy on the frontal or lateral plane.

ASSESS

Visualize the embolized vessels in 3D and their relationship to the other structures with 3DCT HD subtracted.

Distal aneurysm embolization



PLAN

Virtually and dynamically interrogate vessels to understand the relationship with distal aneurysm. Automatically track vessel pathways to anticipate navigation challenges.

GUIDE

Export automatically segmented vessels and landmarks of potential embolization points for augmented-3D fusion guidance. Overlay them on live fluoroscopy on the frontal or lateral plane.

ASSESS

Visualize the embolized vessels in 3D and their relationship to the other structures with High Definition 3DCT.


- 3DCT HD is an option sold separately. Includes 3D XR. Requires AM workstation and Volume Viewer. Available only on Innova IGS 5, Innova IGS 6, Discovery IGS 7 and Discovery IGS 7 OR.
- Split Phase is a GE created customized protocol which requires Volume Viewer Application sold separately.

Needle ASSIST


Improve speed and accuracy

Performing needle procedures in the interventional suite frees up your CT system and provides exceptional access to the patient. However, under fluoroscopic guidance, it may be challenging and time-consuming to find the right entry point and advance the needle while avoiding critical structures.

Needle ASSIST¹⁷ provides real-time visualization of needle positions in the 3D space. It can help medical professionals improve their accuracy, reduce dose and support efficiency efforts when performing needle interventions, while having a limited impact on workflow. The potential reduction in dose and time ultimately helps increase the procedure volume in the angio suite. This ultimately can result in maximizing ROI, free-up CT scan time for improved CT-ROI for diagnostic purposes.




Needle ASSIST enables:



3D SPACE RECONSTRUCTION FOR THE NEEDLE VISUALISATION AT A FULL CBCT WORKFLOW

-98%

RADIATION DOSE¹⁸




TO RECONSTRUCT A NEEDLE IN 3D WITH TWO FLUOROSCOPIC IMAGES WITH AN ACCURACY¹⁹ BETTER THAN

2.5 mm

Stereo 3D

- 1** First fluoroscopic shot and computation of second angulation.
- 2** Second angulation reached thanks to autopositioner and second fluoroscopic shot.
- 3** Validation of registration and automatic needle detection.
- 4** Automatic needle reconstruction and display on the CBCT.

Spine Kyphoplasty



PLAN

With high definition CBCT image the spine with just one fast 3D spin. Determine the optimal pedicle entry point and target defining the needle path directly on oblique CBCT cross sections.


GUIDE

Export the needle trajectories to live fluoroscopy for augmented 3D fusion guidance following C-arm and table movements and use automatically computed bull's eye and progress views for each trajectory. 3D bone rendering helps you visualize mis-registrations and correct for small patient motion from table-side.

ASSESS

Evaluate cement placement with high definition CBCT for immediate quality control of technical success of the procedure.

Pelvic bone osteosynthesis



PLAN

Define optimal needle entry points and pathways directly on oblique CBCT cross sections.

GUIDE

3D fusion guidance all at table-side. Visualize and correct mis-registrations in both translation and rotation.

ASSESS

All at table-side, reconstruct a needle in 3D with 2 fluoroscopic images with accuracy and to review the location of the reconstructed needle on the 3D anatomy.

Optimized for a wide range of specialties

Cardiac and structural heart


Whether it is for percutaneous coronary or valvular interventions, the Innova IGS 6 features transform itself in one of the most efficient monoplane systems. The off-set C-arm allows wide groin coverage without the need to move the L-arm. Combined with the 30 cm detector size it provides open access to the patient at all times & steep angulations, with +/- 55 deg cranial to caudal angulations.



PCI ASSIST^{20,21}
Increased visibility at same dose²²



Quantitative Analysis Package²⁰
Stenosis and Left Ventricle measurements



Calcification Enhancement²³
Improves the visualization of moving contrasted structures



Valve ASSIST²³
Fuse in real time 2D fluoroscopy with 3D anatomy from CT/X-ray

Interventional radiology

The Innova IGS 6 provides positioning flexibility and head to toe coverage necessary for peripheral vascular and abdominal interventions. The 30 cm panel size allows DSA and Cone Beam CT acquisitions on large Field Of View to image the entire anatomy of interest.



InnovaBreeze²⁴
Follow the contrast while looking at subtracted images in real time. Image both legs at the same time with the 30 cm detector.



3DCT HD¹⁰
Allows to image large organs such as liver



Liver ASSIST V.I.²⁵
Simulates injection points in real-time



Vessel ASSIST¹⁴
2D-3D Fusion image guidance for complex CTOs

References

- Innova IGS 6 refers here to IGS 650 configuration
- Optional features, not available on all IGS configurations, refer to your sales representative
- Autoflight features of GEHC's Image Chain requires of GEHC's Interventional X-ray systems, from image acquisition to image processing (IGS 5 and Discovery IGS 7 OP. Not available for sale in all regions.
- Based on competitive analysis, among major players in interventional imaging.
- ASSIST solutions are composed of multiple medical devices. For more information, please refer to GEHC's web site: www.gehealthcare.com/assist. Performance obtained from following publicly available peer reviewed papers: Novel integrated 3DCT and fluoroscopy fusion for LAAC. Value of Image Fusion Coronary angiography for the detection of CABG. Impact of Hybrid rooms with Image fusion on radiation exposure during endovascular Aortic repair. Percutaneous Bone Biopsies: comparison between CBCT and CT guidance. Significant patient radiation exposure reduction during complex liver IR procedures using a new generation angiography imaging room. Comparison of the number of image acquisitions and procedural time required for TACE of Hepatocellular Carcinoma with and without tumor feeder detection SW. In clinical use outcomes will vary depending on the system, settings, clinical task, patient size, anatomical location, clinical practice and ASSIST solutions.
- This type of contrast may not be available in your area. Contact your sales representative for more information.
- Vascular ONWatch study documented in infographic
- GE Internal data
- Statistics for imaging equipment in 2018 in Europe.
- ONWatch is an optional feature, not available in all countries and sold separately.
- mVQ is a customization of image display available on Innova IGS 5, Innova IGS 6, Discovery IGS 7 & Discovery IGS 7 OR. Refer to your sales representative. System is delivered with default settings. Customization requires a GEHC representative.
- For XA modality series, integrated registration currently supports only 3D X-ray Angiography images (stored as CT image storage DICOM objects) acquired with GE mona equipment and reconstructed with the Innova 3DCA application.
- 3DCT HD is an option sold separately. Includes 3DCA. Requires AV workstation and Volume Viewer.
- Optional features on Innova IGS 6
- Option sold separately. Metal artifacts get removed typically by coils and clips in the CBCT FOV.
- Requires AV workstation.
- Vessel ASSIST solution includes Vision 2, VesselIQ Express and Volume Viewer. These applications are sold separately. Not available for sale in all regions.
- Magnetic Resonance Venogram 3D Live Guidance. During Venous Stinus Stenting Treatment and Safety. T Link et al. In J Neurointervention Surg 2017; 10(Suppl 2). Technical note: Bi-plane 3D guidance techniques for neuro-interventional procedures: feasibility and potential benefits. A Santillan et al. In J Neurointervention Surg 2017; 10(Suppl 2). The Statements by GE's customers published here are based on results that were achieved in the customer's unique setting. Results may not be replicable by other GE users.
- Embo ASSIST solution includes FlighPlan for Embolization, Vision 2, VesselIQ Express, Autolone Xpress and requires AV workstation with Volume Viewer, Volume Viewer. These applications are sold separately. Not available for sale in all regions.
- Needle ASSIST solution includes TrackVision 2, stereo 3D and requires AV workstation with Volume Viewer, Volume Viewer. These applications are sold separately. Not available for sale in all regions.
- Based on the dose of the procedure step needed for needle visualization using a CBCT acquisition vs. a Stereo 3D process. Full 3D anatomic information is provided with the CBCT acquisition, while the Stereo 3D process provides specific information for 3D needle visualization. In both cases, the needle visualization is next used to assess its location. The stated dose reduction does not reflect the entire interventional procedure, but rather to a specific step in the procedure. The dose for the CBCT acquisition is from typical exposure settings (Innova CT 407/5, 50ppm, IQ Standard, Normal, Nominal FOV).
- The dose from the Stereo 3D process is from three spatially separated, 2-second fluoroscopic acquisitions, with typical exposure settings (375 pp, IQ Standard, Normal, Max Dose Reduction, Nominal FOV). The dose data for all GE 60007 2-4 seconds, provided for the clinical practice the use of Stereo 3D in clinical practice the radiation dose depending on the clinical task, patient size, anatomical location and clinical practice.
- The accuracy is defined to be the perpendicular distance between the needle tip in the Stereo 3D image and the shaft of the needle in the CBCT image. This accuracy does not reflect the error in the direction parallel to the needle shaft. The perpendicular accuracy was determined by engineering analysis using rigid phantom data. This idealized accuracy of the Stereo 3D reconstruction is obtained with the 2 fluoroscopic images taken at optimal angulation and without table motion at any step of the reconstruction procedure.
- The accuracy is defined to be the perpendicular distance between the needle tip in the Stereo 3D image and the shaft of the needle in the CBCT image. This accuracy does not reflect the error in the direction parallel to the needle shaft. The perpendicular accuracy was determined by engineering analysis using rigid phantom data. This idealized accuracy of the Stereo 3D reconstruction is obtained with the 2 fluoroscopic images taken at optimal angulation and without table motion at any step of the reconstruction procedure.
- Improvement vs. same test without PCI ASSIST. IQ & visibility improvement is measured on Innova IGS530 with thickness using various Plexiglas Phantoms; acquisition parameters and the NEMA spoke wheel tool (ref 11), calculating the ratio of the contrast of the moving wires to the background noise level. The amount of IQ improvement related to PCI ASSIST depends on the acquisition parameter, clinical task, patient size, amount of motion in the image, anatomical location, and clinical practice. Ref 1: A new tool for benchmarking cardiovascular fluoroscopes; S. Balter. Radiation Protection Dosimetry, Vol. 94, No. 1-2 pp. 161-166 (2001).
- Improvement vs. same test without PCI ASSIST. IQ & visibility improvement is measured on Innova IGS530 with Phantoms using various Plexiglas thickness; acquisition parameters and the NEMA spoke wheel tool (ref 11), calculating the ratio of the contrast of the moving wires to the background noise level. The amount of IQ improvement related to PCI ASSIST depends on the acquisition parameter, clinical task, patient size, amount of motion in the image, anatomical location, and clinical practice. Ref 1: A new tool for benchmarking cardiovascular fluoroscopes; S. Balter. Radiation Protection Dosimetry, Vol. 94, No. 1-2 pp. 161-166 (2001).
- PCI ASSIST solution includes SternViz and SternViz2, features of Interventional X-ray systems Innova IGS 5, Innova IGS 6, Discovery IGS 7 and Discovery IGS 7 OR.
- Option sold separately.
- PCI ASSIST solution includes SternViz and SternViz2, features of Interventional X-ray systems Innova IGS 5, Innova IGS 6, Discovery IGS 7 and Discovery IGS 7 OR.
- Improvement vs. same test without PCI ASSIST. IQ & visibility improvement is measured on Innova IGS530 with Phantoms using various Plexiglas thickness; acquisition parameters and the NEMA spoke wheel tool (ref 11), calculating the ratio of the contrast of the moving wires to the background noise level. The amount of IQ improvement related to PCI ASSIST depends on the acquisition parameter, clinical task, patient size, amount of motion in the image, anatomical location, and clinical practice. Ref 1: A new tool for benchmarking cardiovascular fluoroscopes; S. Balter. Radiation Protection Dosimetry, Vol. 94, No. 1-2 pp. 161-166 (2001).
- Improvement vs. same test without PCI ASSIST. IQ & visibility improvement is measured on Innova IGS530 with Phantoms using various Plexiglas thickness; acquisition parameters and the NEMA spoke wheel tool (ref 11), calculating the ratio of the contrast of the moving wires to the background noise level. The amount of IQ improvement related to PCI ASSIST depends on the acquisition parameter, clinical task, patient size, amount of motion in the image, anatomical location, and clinical practice. Ref 1: A new tool for benchmarking cardiovascular fluoroscopes; S. Balter. Radiation Protection Dosimetry, Vol. 94, No. 1-2 pp. 161-166 (2001).
- Results obtained at Valley Baptist Medical Center, Hartmann, TX, USA for stroke patients receiving mechanical thrombectomy and enrolled in the STROKE registry. The statements by GE's customers described here are based on results that were achieved on the customer's unique setting. These results are specific to the patient size, case mix, there on. We do not guarantee that other customers will achieve the same results.
- Endovascular Treatment Outcomes Using the Stroke Triage Education, Procedure Standardization and Technology (STEPS-1) Program. Hassan et al. In Interv Neuro radiol. 2018 Feb;24(1):51-56. doi: 10.1177/1591019917740100.
- EXTEND-IA Trial. Campbell BC et al: Endovascular therapy for ischemic stroke with perfusion-imaging selection. N Engl J Med 2015; 372:1009-1018. DOI:10.1056/NEJMoa1414792
Endovascular treatment: improved the functional outcome at 90 days, with more patients achieving functional independence (score of 0 to 2 on the modified Rankin scale, 71% vs. 40%, p=0.01).
- Requires an AV workstation with Volume Viewer and Volume Viewer Innova. These applications are sold separately.



GE Healthcare

About GE Healthcare

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE:GE) works on things that matter – great people and technologies taking on tough challenges. From medical imaging, software & IT, patient monitoring

GE Healthcare is a leading global medical technology and digital solutions innovator. GE Healthcare enables clinicians to make faster, more informed decisions through intelligent devices, diagnostic pharmaceuticals, data analytics, applications and services, supported by its Edison intelligence platform. With over 100 years of healthcare industry experience and around 50,000 employees globally, the company operates at the center of an ecosystem working toward precision health, digitizing healthcare, helping drive productivity and improve outcomes for patients, providers, health systems and researchers around the world.

Follow us on Facebook, LinkedIn, Twitter and Insights or visit our website www.gehealthcare.com for more information

GE Healthcare
Chalfont St.Giles,
Buckinghamshire,
UK



GE Healthcare



@GE Healthcare



GE Healthcare



GE Healthcare

Data subject to change.

Marketing Communications GE Medical Systems

Société en Commandite Simple au capital de 65.146.245 Euros

283 rue de la Minière – 78530 Buc Cedex France

RCS Versailles B 315 013 359

A General Electric Company, doing business as GE Healthcare

GE, GE Monogram, Innova IGS 6, AutoRight and Innova Breeze are trademarks of General Electric Company.

Innova IGS 6 and products mentioned in this material cannot be marketed in countries where market authorization is required and not yet obtained. Refer to your sales representative.

GE Healthcare, Europe
Headquarters Buc, France
+33 800 90 87 19

GE Healthcare, Asia Pacific
Tokyo, Japan
+81 42 585 5111

GE Healthcare, Middle East and Africa
Istanbul, Turkey
+90 212 36 62 900

GE Healthcare, ASEAN
Singapore
+65 6291 8528

GE Healthcare, North America
Milwaukee, USA
+1 866 281 7545

GE Healthcare, China
Beijing, China
+86 800 810 8188

GE Healthcare, Latin America
Sao Paulo, Brazil
+55 800 122 345

GE Healthcare, India
Bangalore, India
+91 800 209 9003

Imagination at work

©2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown hereby, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE and the GE Monogram are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

JB83509XX

gehealthcare.com